



# Pipetting Proficiency Test Method




# Background

- Pipetting variance can lead to errors which can:
  - Result in unrecognized wrong results
  - Failure to pass a good product (as in QC)
  - Failure to fail a poor product (as in QC)
- Pipetting variance can be attributed to:
  - Poorly seated tips
  - Variable viscosities not taken into account
  - Profusion (residual volume on exterior of tip) carry over
  - Pipette out of calibration (dropped or overly used)
  - Inexperience or need for retraining



# Definition

Despite training or prior knowledge, it is possible to drift from uniform practice (as in a golf swing). While pipetting may be effective, it may not yield an accurate and/or consistent result.





# Purpose

- This Test Method is designed to characterize pipetting efficiency across the organization.
- Depending on results, there may or may not be a need for training/re-training.
- Uniform pipetting practices are essential to the success of Salimetrics and to our customers.
- This test method, or a revision, is intended to be developed for our customers to help them ensure that their pipetting methods are correct.



# Approach

- Who: all individuals who routinely operate pipettes are being asked to participate in this evaluation.
- What: Fluids of varying viscosity will be pipetted by each operator into microwell plates.
- When: within the next 3 days. The procedure should only take a half hour.
- Where: at your own lab bench, with your current pipette
- Why – ibid



# Methods

- Fluids of varying viscosities (xanthan gum, water, EtOH) containing yellow dye were developed for this study
- The fluids will be pipetted in multiple replicates by each operator into microwell plates
- All additions are first pass only – no repeats

## **Single channel pipetting**

10 uL (viscous) into dry wells, add 100 uL water, 48 replicates

50 uL (viscous) into dry wells, add 50 uL water, 48 replicates

100 uL (viscous) into dry wells, 48 replicates

## **Multi channel pipetting**

50 uL (viscous) into dry wells, add 50 uL water, 48 replicates

100 uL (viscous) into dry wells, 48 replicates

200 uL (25% EtOH) into dry wells, 48 replicates



# Materials

Each operator will be provided:

4 dye fluids made with McCormick's Yellow Dye

- 10 uL fluid (1 mL), dilute dye 58-fold in xanthan gum solution
- 50 uL fluid (8 mL), dilute dye 291-fold in xanthan gum solution
- 100 uL fluid (15 mL), dilute dye 580-fold in xanthan gum solution
- 200 uL fluid (15 mL), dilute dye 1160-fold in 25% EtOH solution

6 uncoated microwell plates

## **Xanthan Gum Solution**

0.95 g/L Xanthan gum in water

0.35 g/L methyl p-hydroxy benzoate (preservative)

## **25% EtOH Solution**

25% EtOH in water


0.35 g/L methyl p-hydroxy benzoate (preservative)



# Pipetting Patterns – Single channel


## 10sc (single channel)

Add 10 uL dye fluid to each of well dry wells.  
Add 100 uL water to each column of wells with multichannel pipetter.  
Read absorbance at 450 – 492. Export to F: drive




## 50sc (single channel)

Add 50 uL dye fluid to each of well dry wells.  
Add 50 uL water to each column of wells with multichannel pipetter. Read absorbance at 450 – 492. Export to F: drive



## 100sc (single channel)

Add 100 uL dye fluid to each of well dry wells.  
Read absorbance at 450 – 492. Export to F: drive








# Pipetting Patterns - Multichannel


## 50mc (multi channel)

Add 50 uL dye fluid to each of well dry wells with all eight channels, left to right. Add 50 uL water to each column of wells with multichannel pipetter. Read absorbance at 450 – 492. Export to F: drive




## 100mc (multi channel)

Add 100 uL dye fluid to each of well dry wells with all eight channels, left to right. Read absorbance at 450 – 492. Export to F: drive



## 200mc (multi channel)

Add 200 uL dye fluid to each of well dry wells with all eight channels, left to right. Read absorbance at 450 – 492. Export to F: drive





# Analysis (GAGE R&R)

The results will be uploaded into MiniTab and analyzed by ANOVA for:

- well-to-well (repeatability)
- operator-to-operator (reproducibility)

**Gage R&R Study (Crossed)**

C1	StdOrder	Part numbers:	Parts	Gage Info...
C2	RunOrder	Operators:	Operators	Options...
C5	Food Temperature	Measurement data:	Food Temperature	Conf Int...
		Method of Analysis:	<input checked="" type="radio"/> ANOVA <input type="radio"/> Xbar and R	Storage...

Buttons: Select, Help, OK, Cancel

Callouts:

- 10sc, 50sc, 100sc, 50mc, 100mc, 200mc (pointing to Part numbers)
- OD's (pointing to Operators)
- PP1 - PP20 (pointing to Method of Analysis)

The manual and automatic pipetter results may need to be handled separately.

# Pipetting Proficiency Evaluation


## Methods

- Artificial saliva fluids with yellow dye were developed to deliver an OD of 1.0 after pipetting 10, 50, 100 or 200 uL into microwell plates.
- Plates were read on a plate reader and the data uploaded to the f: drive, keeping the identity of the operators anonymous.
- Data was analyzed, by MiniTab: ANOVA, Box/Whisker, Dotplot and GAGE R&R

## Example – single channel

### 10sc (single channel)

Add 10 uL dye fluid to each of well dry wells.  
Add 100 uL water to each column of wells with multichannel pipetter.  
Read absorbance at 450 – 492. Export to F: drive




# Pipetting Proficiency Evaluation

## Example – multi channel

### 100mc (multi channel)

Add 100 uL dye fluid to each of well dry wells with all eight channels, left to right.

Read absorbance at 450 – 492. Export to F: drive



## Results

Thirteen (13) operators each pipetted dye solutions into six plates (48 wells each) for a total of 288 wells

- 10 uL, 50 uL, 100 uL single channel
- 50 uL, 100 uL, 200 uL multichannel

Participants were from QC, Testing, Manufacturing and R&D

# 10 uL, single channel, worst case fluid, viscous, high dye concentration, small volume






<u>Within Operator</u>	
Avg CV =	8.2%
Stdev =	2.6%

<u>Operator - Operator</u>	
CV =	8%

<u>Pooled (Total)</u>	
CV =	11%

<u>Nested ANOVA:</u> 10 uL SC –vs. Operator, Replicate Variance Components	
Source	
Operator-Operator	46 %
Well-Well	54 %



Example Fluids: Conj – E2, E3, Prog, Testo  
Sample – α Amylase, Cort (25), Testo (25)

# 50 uL, single channel, viscous

## Within Operator

Avg CV = 8.3%

Stdev = 2.0%

## Operator - Operator

CV = 9%

## Pooled (Total)

CV = 12%


## Nested ANOVA:

10 uL SC –vs. Operator, Replicate  
Variance Components

Source


Operator-Operator 53 %

Well-Well 47 %



Example Fluids : Conj – BC, E1

Sample – Prog, DHEA



# 100 uL, viscous, single channel

## Within Operator

Avg CV = 2.0%

Stdev = 1.0%

## Operator - Operator

CV = 4%

## Pooled (Total)

CV = 4%


## Nested ANOVA:

10 uL SC –vs. Operator, Replicate  
Variance Components


Source

Operator-Operator 75 %

Well-Well 25 %



Example Fluids : Sample – E1, E2, E3, DHEA-S



# 50 uL, viscous, multi channel

## Within Operator

Avg CV = 7.6%

Stdev = 2.2%

## Operator - Operator

CV = 9%

## Pooled (Total)

CV = 12%



## Nested ANOVA:

10 uL SC –vs. Operator, Replicate  
Variance Components

Source

Operator-Operator 58 %

Well-Well 42 %





# 100 uL, viscous, multi channel

## Within Operator

Avg CV = 1.7%

Stdev = 0.7%

## Operator - Operator

CV = 1.7%

## Pooled (Total)

CV = 2.5%



## Nested ANOVA:

10 uL SC –vs. Operator, Replicate  
Variance Components

Source

Operator-Operator 44 %

Well-Well 56 %



Example Fluids :

Sample – E1, E2, E3 for HX10 testing

# 200 $\mu$ L, low viscosity, multi channel



## Within Operator

Avg CV = 1.1%  
Stdev = 0.3%

## Operator - Operator

CV = 0.5%

## Pooled (Total)


CV = 1.3%

## Nested ANOVA:


10  $\mu$ L SC –vs. Operator, Replicate  
Variance Components

Source

Operator-Operator 11 %  
Well-Well 89 %



Example Fluid: TMB substrate





# Summary and Conclusions

- Levels of error are within expected ranges for manual pipetting.
- Larger volumes typically yield better precision. This is true in automated systems also.
- Pipette types (Hamilton, Finnpipette, Finnpipette 2, BioHit) were not factors (data not shown)
- Testing method was easy to implement, follow and execute.
- Remind operators to ensure tips are well seated and to minimize profusion (already common practice).